


MISE A JOUR LE :	PAR :	NATURE DES MODIFICATIONS :
29/12/2021	G. DELVALLEZ	English translation
RENOUVELLEMENT LE :	PAR :	

1. SAMPLING CONDITIONS

- Morning urine or continence of at least 3 hours
- Sterile container
- Desinfectant (Dakin)
- Hand washing soap



2. SAMPLE



- Hand washing with soap
- Male: Dakin disinfection of the meatus and glans orifice
- Woman: Dakin disinfection of the opening of the urethral and vaginal meatus
- Elimination of the first drops of urine
- Collection of the mid-stream in the sterile container
- Minimum quantity: 20 ml

3. TRANSPORTATION

- The sample should be sent to the laboratory within 2 hours.
Otherwise, it can be stored in the refrigerator (2-8 ° C) for 6-8 hours.
- Indicate the treatment prescribed by the doctor (antibiotics)

**La mise à jour de ce document est garantie sur KaliWeb.
Veuillez régulièrement à l'actualisation de vos éditions papier.**